




# Development of custom, low-cost multiscale microscopic imaging systems for neuroscience applications

Ashly Jose\*<sup>1</sup>, Zahra Laouby<sup>1,3</sup>, Stanley Tang<sup>2</sup>, Yukti Vyas<sup>1,3</sup>, Marco Bonesi<sup>1</sup>, Claude Aguergaray<sup>2</sup>, Neil Broderick<sup>1</sup>, Frederique Vanholsbeeck<sup>1</sup>, Juliette Cheyne<sup>1,3</sup>

<sup>1</sup>The Dodd-Walls Centre for Photonic and Quantum Technologies, Department of Physics, The University of Auckland, New Zealand. <sup>2</sup>Photon Factory, The University of Auckland, New Zealand, <sup>3</sup>Physiology Department and Centre for Brain Research, The University of Auckland, New Zealand. \* ajos486@aucklanduni.ac.nz



## Introduction

- It is now possible to image brain cell activity in awake behaving animals; however, commercial systems are expensive, prohibiting their widespread use.
- To overcome this, we are developing cost-effective and easily customizable imaging tools (mesoscope and two-photon microscope) to visualize structural and functional details of the brain.
- Mesoscope - fluorescence imaging for studying meso-scale activity from diverse regions of the brain for understanding the functional properties and correlations of widely separated neurons in the brain.
- Two-photon microscope (TPM) - high-resolution imaging using ultra short pulses from a passively mode-locked fibre laser for studying micro-scale activity from individual neurons.

## Theory and Design

- One-photon excitation: a single photon excites the fluorophore.
- Two-photon excitation: two low-energy photons combine to excite the fluorophore.

- The average rate of two photon excitation per molecule :

$$\text{Rate} \equiv \frac{1}{2} \delta \langle I^2 \rangle$$

where,  $\delta$  is two photon absorption cross section and  $\langle I^2 \rangle$  is the time average of second order local laser intensity




Fig 1. Principle of single-photon and two-photon excitation [3].

### Why Mesoscope?

Large Field of View  
Image entire mouse cortex  
Examine multiple cortical areas simultaneously

### Why TPM?

Localized excitation/reduced phototoxicity  
Good depth penetration  
Excellent resolution (subcellular)

## Current Directions

- Optimization of the mesoscope for *in vivo* fluorescence imaging in the entire cortical region of the mouse brain.
- Calibrating the TPM for *in vivo* fluorescence imaging in the auditory cortex of mouse brain
- Imaging response to sound in the auditory cortex of ASD mice compared to wild type.
- Genetically encoded  $\text{Ca}^{2+}$  indicators act as neuronal activity reporters - GCaMP for mesoscope and r-GECO for TPM.

	TPM	Mesoscope
Excitation	Mode locked pulsed laser	LED source
	Wavelength - 1030nm	Wavelength - 470nm
FOV	<1mm x 1mm	>10mm x 10mm
Detection	Photo multiplier tube	CMOS camera




Fig 3. Images acquired using mesoscope (a) *In vitro* calcium imaging of brain slices (b) *In vivo* calcium imaging of neonate rat brain.

## Conclusion

- We have built a fully customized TPM with 3D movement and rotation, that has given promising results for *in vitro* calcium imaging of brain slices.
- We have developed a low-cost, reverse tandem lens mesoscope for imaging cortical region of the mouse brain.
- These techniques will be combined to better understand neuronal circuit changes in rodent models of disease, including Autism Spectrum Disorder (ASD) and Spinal Cord Injury.

## References

1. Benninger, R. K., & Piston, D. W. (2013). Two-photon excitation microscopy for the study of living cells and tissues. *Current protocols in cell biology*, 59(1), 4-11.
2. Rubart, M. (2004). Two-photon microscopy of cells and tissue. *Circulation research*, 95(12), 1154-1166.
3. P. Bowen, H. S., A. Runge, R. Provo, and N. G. R. Broderick. Mode-locked femtosecond all-normal all PM Yb-doped fiber laser at 1060nm. *Optics Communications* 364, 181-184 (2017).

## Acknowledgements

I gratefully acknowledge the Dodd-Walls Centre for Photonic and Quantum Technologies and Faculty Research Development Fund for supporting this research.




Fig 2. The 3D SketchUp designs of the experimental setup for a) the Mesoscope and b) the Two-Photon Microscope (with photo inset).